



**WALL CABINET
FOR AUTOMATED
EXTERNAL
DEFIBRILLATOR
(AED)**



Aivia® 'S'

OVERVIEW

The AIVIA® «S» is designed for an indoor positioning. Combining strength and lightness, the AIVIA® «S» is composed of 2 distinct parts:

- The polycarbonate cover,
- The ABS chassis.

These materials are well known in the industry, giving excellent shock resistance properties.

Access to the AED is protected by security seals.



LABELING

The survival sequence and three clear pictograms with instruction for use are shown on the front face.



18501 Pines Blvd, 3rd floor
Pembroke Pines, FL 33029, USA