

Certificate of Analysis


CJC-1295 without DAC 2 mg

also known as Modified GRF (1-29)

Compound : CJC-1295 no DAC **Client** : SW Peptides
Lot number : SW-A004-02 www.swpeptides.com
Analysis date : 2022-12-23
Purity % : 99.9%
Method : Mass Spectrometry and UV

PubChem CID : 91976842

<https://pubchem.ncbi.nlm.nih.gov/compound/91976842>



Time (min)	Area	%Area	
16.60	1.86E+01	0.4	
16.99	4.48E+03	99.6	CJC-1295 no DAC

Peak identity confirmed by mass spectrum evaluation
 Purity confirmed using UV detection
 Expected mass : 3365.89 g
 Measured mass : 3365.89 g
 Molecular weight confirmed

Note: Injectable peptides contain salts and sugars to aid in solubility and act as pH buffers. These are not normally detected using UV and are not considered impurities.

Analysis Performed by
 Ken Pendarvis, ChE
 Analytical Chemist
 MZ Biolabs
contact@mzbiolabs.com



2022-12-23